

Release Information

Release Date: December 2024

Release Build: 1.19.7-6 (patch)

GALE App Fixes & Features (Overview)

Main Features

- Customizable Intake and Wrap-up forms
- Customizable forms to send and gather BA data or health data
 - Option to send to external analytic applications
 - Option to send to email
- Restructure Patient Profiles for Medical Record support
- Federated Authentication login using a Care Provider account
- Update Zoom SDK and library with a fix to support external cameras

Major Bug Fixes

- Use external cameras while on a Zoom call
- Omron Blood Pressure Cuff
- Spirometer Crash
- QT Medical Crash
- Kosmos stability
- Patch: Missing Clinics

Web App Fixes & Features (Overview)

Deployment Management System (DMS)

- Federated Authentication login
- Customizable Intake and Wrap-up forms

Insights / Call Center

- Add additional profile details to the Pre-visit popup

- Allow Providers with access to multiple Care Providers to see calls from all Clinics they are subscribed to regardless of which Call Center they are logged in to

NEW: Remote Patient Management Portal (XMR)

- FHIR-based patient management portal
- Manage Patient
 - Encounters
 - Observations
 - Images and Documents
- Administrative forms for BA

Wellness API

Known Issues

- DMS Forms:
 - Cannot change the order of pages or questions
- Measurements are not displayed in the Zoom sensor panel if Insights is accessed while the sensor is running. Restarting the panel will resolve the issue.
- Potential for providers to be asked for a passcode while joining Zoom. Refresh the Call Center to rectify
- Long insights load times reported
- Green Zoom screen on start
- Cannot send Visit Summary via MMS/WhatsApp through GALE

Crashes

- Popup fragment crashes that do not affect the user experience

Resolved Issues

All Features:

All Bugs: