

A basic set of tests for validating the Athena-GALE solution.

Test 1: Basic Workflow with Video Visit

App	Description	Result	Notes
	Timestamp of Encounter Start: Tablet Name/ID:		
GALE	Log into Athena (specify)	n/a ▾	
GALE	Search for a patient	n/a ▾	
GALE	Load and check in Patient (specify)	n/a ▾	
GALE	Select "Visit Reason" (specify)	n/a ▾	
Athena	Validate encounter is created	n/a ▾	
GALE	Take vitals (specify)	n/a ▾	
GALE	Take images (specify)	n/a ▾	
GALE	Place a call. Fill out intake if applicable (specify clinic/practice)	n/a ▾	
CC	Answer call from Call Center	n/a ▾	
CC	Validate intake/measurements in Zoom chat	n/a ▾	
Insights	Open Insights <i>Validate intake, visit reason, vitals</i>	n/a ▾	
GALE	Take a measurement during a call (specify). <i>Watch for vitals and messages coming through Zoom chat</i>	n/a ▾	
Insights	Add and save a clinical note <i>(if joined from Zoom, this option is disabled)</i>	n/a ▾	
GALE	End Call	n/a ▾	
GALE	End Encounter or log out. <i>Note toast messages/popups</i>	n/a ▾	
	Validate Encounter	n/a ▾	
Athena	Validate encounter on Athena	n/a ▾	

Test 2: Basic Workflow without Video Visit

App	Description	Result	Notes
	Timestamp of Encounter Start: Tablet Name/ID:		
GALE	Log into Athena (specify)	n/a ▾	
GALE	Load and check in Patient (specify)	n/a ▾	
GALE	Select "Visit Reason" (specify)	n/a ▾	
Athena	Validate encounter is created	n/a ▾	
GALE	Take vitals (specify)	n/a ▾	
GALE	Take images (specify)	n/a ▾	
GALE	End Encounter or log out. <i>Note toast messages/popups</i>	n/a ▾	
	Validate Encounter	n/a ▾	
Athena	Validate encounter on Athena	n/a ▾	

Test 3: Guest Encounter - Video Visit workflow without logging in

App	Description	Result	Notes
	Timestamp of Encounter Start: Tablet Name/ID:		
GALE	Dismiss Login dialogue	n/a ▾	
GALE	Take vitals (specify)	n/a ▾	
GALE	Take images (specify)	n/a ▾	
GALE	Place a call. Fill out intake if applicable (specify clinic/practice)	n/a ▾	
CC	Answer call from Call Center	n/a ▾	
CC	Validate intake/measurements in Zoom chat	n/a ▾	
Insights	Open Insights <i>Validate intake, visit reason, vitals</i>	n/a ▾	
GALE	Take a measurement during a call (specify). <i>Watch for vitals and messages coming through Zoom chat</i>	n/a ▾	
Insights	Add and save a clinical note <i>(if joined from Zoom, this option is disabled)</i>	n/a ▾	
Insights	Share and download summary. <i>Validate notifications</i>	n/a ▾	
GALE	End Call	n/a ▾	
GALE	End Session <i>Note toast messages/popups</i>	n/a ▾	

Test 4: Individual measurements

Validate each encounter on Athena after each test. This set of tests is to ensure that each individual measurement is configured correctly. Note any unexpected results and contact 19Labs with specific details.

Vitals & Images *with* a video visit

App	Description	Result	Notes
	Timestamp of Encounter Start: Tablet Name/ID:		
GALE	Log in / load patient / take only otoscope images / place call / clinical note + save / end encounter	n/a ▾	
GALE	Log in / load patient / take only blood pressure / place call / clinical note + save / end encounter	n/a ▾	
GALE	Log in / load patient / take only pulse oximeter / place call / clinical note + save / end encounter	n/a ▾	
GALE	Log in / load patient / take only temperature / place call / clinical note + save / end encounter	n/a ▾	

Vitals & Images *without* a video visit

App	Description	Result	Notes
	Timestamp of Encounter Start: Tablet Name/ID:		
GALE	Log in / load patient / take only otoscope images / end encounter	n/a ▾	
GALE	Log in / load patient / take only blood pressure / end encounter	n/a ▾	
GALE	Log in / load patient / take only pulse oximeter / end encounter	n/a ▾	
GALE	Log in / load patient / take only temperature / end encounter	n/a ▾	

Add as many additional tests as you can think of. Give them a title, and outline the steps you took and the result.

Test 5: Additional Tests

[illegible]